



Post Infusion Cooling Times

Chemotherapy Regimen	Post Infusion Cooling Time (minutes)
Paclitaxel 175mg/m ² (q2 week)	120
Docetaxel/Cyclophosphamide	120
Docetaxel/Carboplatin	120
Paclitaxel (weekly dose)	90
Doxorubicin/Cyclophosphamide	180 (European clinical study recommendations)
Docetaxel/Carboplatin/Herceptin	120 post Carboplatin
Docetaxel/Carboplatin/Herceptin/Perjeta	120 post Carboplatin
Docetaxel/Cyclophosphamide/Herceptin	120 post Cyclophosphamide
Abraxane	90
Abraxane/Gemzar	120
Carboplatin or Cisplatin/Gemzar	120
Cytosan/Methotrexate/5-FU	120
Eribulin	180
Paclitaxel 175mg/m ² + Carboplatin (q3 week)	120
Docetaxel mono-therapy (q3 week)	120
Other Taxane-based regimens (not weekly)	120
Other Anthracycline-based regimens	180

PRO-20170801-03-US 09/2018

DigniCap Clinical Support Center

www.dignicap.com/helpdesk | 877-350-2150 | support@dignicap.com