

TECHNICAL SPECIFICATIONS

Dimensions	23" wide x 25" deep x 43" tall including wheels
Weight	Approx. 165 pounds with empty tank
Mobility	Machine is on wheels and can be easily moved
Power Supply	115 VAC
Fuse 115 VAC	16 A
Alarms	Safety system Low coolant level Temperature out of range DigniCap not connected Pause (when pause > 8 min.) Desired temperature in tank not reached Tank temperature too high Tank temperature too low (below -10° C (14° F)) Malfunctioning sensor Malfunctioning sensor cable contact Service indicator System failure
Operating Conditions	Temperature: 64 - 86° F Humidity range: 30 - 90% RH Altitude: ≤ 6561 feet
Refrigeration Unit	Fully hermetically sealed unit, using CFC-free R404A refrigerant
DigniCool (Coolant)	Diluted monopropylene glycol Low toxicity under normal conditions of handling and use Combustible material and a slight eye irritant MSDS Available
Coolant Tank	2.25 Gallons
Cap Material	Silicone inner cap Neoprene outer cap
Sensors in Cap	2 temperature sensors (front & back) 1 safety sensor